



AliMed ONE® Wireless Motion Sensor

Instructions for Use, Setup, and Care

Intended Use

The AliMed ONE Wireless Motion Sensor pairs with the Alarm Base to provide caregiver notification of patient falls or elopement as part of an overall fall management program. The Motion Sensor can be mounted to a wall for wandering protection anywhere in the room.

Key Features

- 1. Universal Status Indicator**
Illuminates at varying times: when Motion Sensor is turned on, when pairing with Alarm Base is complete, or when infrared beam is broken
- 2. Infrared Motion Sensor**
Infrared beam detects motion
- 3. On/Off Switch** Turns Motion Sensor on or off
- 4. Pair Status Indicator**
Illuminates when Alarm Base is in range for pairing
- 5. Battery Compartment Button**
Press to open battery compartment
- 6. Universal Wall Mount and Bracket** Holds Motion Sensor and mounts to wall




Figure 1



Figure 2




Figure 3

Setup and Operating Instructions

The Motion Sensor may be set up by Nurse or Biomedical staff.

- Motion Sensor requires one (1) 9V battery to operate (not included).
- Press the Battery Compartment Button (Fig. 3) to open the battery compartment.
- Align the male and female connectors on the top of the battery to the connectors attached to the unit. Place in battery compartment.
- The Motion Sensor is now ready for pairing.
- Close the battery compartment door until it snaps into place.

Pairing the Motion Sensor with the Alarm Base:

- Ensure the On/Off Switch is set to the "OFF" position (Fig. 2).
- Press Pair Button on the Alarm Base.
- The Base will beep and the first open channel will blink green and red indicating channel is ready for pairing.
- Place Motion Sensor next to the Base matching PAIR on the Motion Sensor to PAIR on the Base (Fig. 4). *Note: Pair Status Indicator on Motion Sensor (Fig. 3) will illuminate when in range of the Base.*
- If pairing is successful, the Base will beep and the assigned channel will blink green once. The Universal Status Indicator (Fig. 1) on Motion Sensor will blink three (3) times.
- The Motion Sensor is now paired with the Base. The assigned channel on the Base will continuously blink at approximately 5-second intervals indicating Motion Sensor is in use.




Match PAIR on Motion Sensor to PAIR on Base

Figure 4

Wall Mounting Instructions

1. Secure the bracket to the wall with the screws provided (Fig. 3).
2. Align mounting slot on Motion Sensor with mounted bracket.
3. Slide Motion Sensor from the left until secured in place on the mounted bracket (Fig. 5).
4. Adjust Motion Sensor angle as needed.
5. The Motion Sensor is now ready to use.



Slide Motion Sensor from the left until secured in place on bracket

Figure 5

Care

1. Wipe with soft, non-abrasive cloth using a 1:10 bleach solution (contains 0.5% chlorine concentration) or follow your facility's cleaning protocol.
2. Do not immerse in fluids.
3. After cleaning, check unit is functioning properly.

Warning

Failure to comply with warning may result in injury or death. This device is not suitable for all individuals. Other devices may be required. This device is not a substitute for visual monitoring by caregiver. The manufacturer does not claim that this device will stop elopement and/or stop falls. It is a supplemental device to be used as part of a facility's overall fall management program. This product is to be used by a licensed caregiver only. As such, it is the entire responsibility of the caregiver to ensure that the system is properly installed and tested. Further, the system is not designed to replace good caregiving practices including, but not limited to:

- a.) Direct patient supervision
- b.) Adequate training for personnel for fall prevention and elopement
- c.) Testing of the system before each use

Report any serious incidents to AliMed and the appropriate regulatory body.

Limited Warranty: Seller warrants that the goods sold will be free from defects of workmanship and manufacturing for a period of one year from the date of sale.

Item Numbers

REF #713817



Manufactured for:



AliMed® Inc.

297 High St., Dedham, MA 02026 USA
PH 800.225.2610 • FAX 800.437.2966



MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany



Environmental Conditions

	Operational	Storage	Transport
Temperature	10°-40°C 50°-104°F	-25°-70°C -13°-158°F	-25°-70°C -13°-158°F
Humidity	30%-75% Non-condensing	0%-95% Non-condensing	0%-95% Non-condensing
Pressure	700-1060 hPa	700-1060 hPa	700-1060 hPa

CONFORMS TO AAMI STD ES60601-1, IEC STD 60601-1-6
CERTIFIED TO CSA STD C22.2 NO. 60601-1, NO. 60601-1-6

Symbol Glossary

	Manufacturer
	Batch Code/Lot Number
	Catalog Number/Item Number
	Serial Number
	European Conformity Mark
	Authorized representative in the European Community
	Tested/certification for compliance to published industry standards
	Consult instructions for use
	Caution: Read all warnings and precautions in instructions for use
	DO NOT THROW IN TRASH
	Indicates FCC has certified device
	Indicates compliance to RoHS and is free of substances banned by RoHS
	Medical Device

AliMed.com

#7721REV0321 CCO1067

©2021 AliMed® Inc. All rights reserved.

Made in China