

Establish a culture of routine for safer patient handling

Follow the Three "A's"

- ✓ **Assess** your department needs, your staffing levels, and your patient population for mobility level and assistance required
- ✓ **Acquire** the right transfer devices for your specific point-of-care needs
- ✓ **Access** to transfer devices makes it easy for staff compliance

Key Takeaway

Always use some type of transfer tool—something is better than nothing

Choose the right mix of devices for safer handling

Lateral Transfer Aids

Device Types: Air-Assisted • Rollboards • Shifters • Transfer Sheets

Department Use: Emergency Room • Operating Room • Imaging Suite • Patient Room

Mobility Level: No to Limited Mobility / Dependent



Limited Mobility Transfer Aids

Device Types: Full-Body Lifts • Sit-to-Stand Lifts • Transfer Boards

Department Use: Operating Room • Patient Room • Bathroom

Mobility Level: No to Limited Mobility / Fully Dependent



Assisted Mobility Transfer Aids

Device Types: Sit-to-Stand Lifts • Transfer Boards • Gait Belts

Department Use: Patient Room • Inter-Departmental Transfer

Mobility Level: Limited to Assisted Mobility / Moderately to Minimally Dependent



Repositioning Aids

Device Types: Turners • Boosters • Air-Assisted

Department Use: Patient Room

Mobility Level: No to Limited Mobility / Fully Dependent



Lateral Patient Transfer: Pull Force 101



Patient Weight Significantly affects the amount of friction created and, in turn, the pull force required—**weight must always be considered during transfers**

Friction The force that acts to resist sliding between two surfaces that are touching, which can slow down or stop the motion of the object and **increase as weight is added**

Don't underestimate RISK

Remember—more injuries occur when transferring **average size patients!**

Pull Force The amount of force needed to move a patient laterally—the **greater the pull force required, the greater the risk for staff injury.**

Lessen injury risk: know what patient handling tools can handle

Air-Assisted



Lowest Pull Force/
Lowest Injury Risk

Maximum stability and comfort

- Best overall option for manual lateral transfer
- Cushion of air significantly reduces friction, which drastically reduces the pull force required compared to other modalities
- Requires only 2 staff to assist
- For any size patient

Rollboards



Lower Pull Force/
Lower Injury Risk

Safely bridges uneven gaps

- Lightweight and easy to handle
- Friction-reducing mechanics and materials help reduce push/pull efforts
- Requires minimum of 2-3 staff
- Models available for average and bariatric size patients
- Hangable for easy access and storage

Shifters



Medium Pull Force/
Medium Injury Risk

Transfer Sheets



Higher Pull Force/
Higher Injury Risk

Emergency situations only

- Use only if other devices are unavailable
- Lower friction option than traditional linen draw sheets
- Often requires more than 4 staff
- Models available for average and bariatric size patients
- Stores anywhere for easy access

Linen Draw Sheets



Highest Pull Force/
Highest Injury Risk

Not recommended

- Only use with another assistive device to reduce injury risk

