

Protect patients against pressure and shear



Head Donuts

Available in Adult Standard, Adult High, Adult Wide, Pediatric, Pedi/Neonate, and Neonatal



Body Positioners
Chest Rolls shown



Upper Extremities
Contoured Armboards shown



Lower Extremities
Heel Cups shown



Positioning Sets
Basic Positioning Set shown

■ 100% viscoelastic gel



■ Impervious to fluids



■ Easy to clean



AliGel Positioners

AliGel Positioners remain an industry standard due to the pressure-relieving properties of a unique gel formula that mimics fatty tissue and moves with the patient to reduce shear and help prevent pressure ulcers. Additionally, AliGel absorbs shock, stretches, and retains its shape better than less dense gel positioners. Made of 100% dry viscoelastic polymers.

For all your patient positioning needs:

- Head and Neck
- Body
- Upper Extremity
- Lower Extremity
- Overlays
- Pediatric
- Frame Pads
- Positioning Sets

See p. 2 for complete listing.

AliGel™ Positioners

DESCRIPTION	SIZE	ITEM #
HEAD AND NECK		
DONUT	ADULT	#95-118
DONUT	ADULT, HIGH	#938139
DONUT	ADULT, WIDE	#938140
DONUT	PEDIATRIC	#95-126
DONUT	PEDI/NEONATE	#938141
DONUT	NEONATAL	#95-124
HORSESHOE	ADULT	#95-115
HORSESHOE	ADULT, HIGH	#938142
HORSESHOE	PEDIATRIC	#95-127
HORSESHOE	PEDI/NEONATE	#938143
HORSESHOE	NEONATAL	#95-125
PRONE HEAD PAD	ADULT (L)	#925815
PRONE HEAD PAD	PEDIATRIC (M)	#931677
SUPINE HEAD PAD	ADULT	#925814
LATERAL HEAD PAD		#95-938
CONTOURED HEADREST		#930137
OPHTHALMIC HEADREST	ADULT	#933000
OPHTHALMIC HEADREST	PEDIATRIC	#933001
OPHTHALMIC CRADLE	ADULT	#932448
OPHTHALMIC CRADLE	PEDIATRIC	#932449
FOREHEAD PROTECTOR		#938144
BODY		
CHEST ROLL	2"W x 6"L x 2"H	#921799
CHEST ROLL	2"W x 8"L x 2"H	#925820
CHEST ROLL	2"W x 10"L x 2"H	#925821
CHEST ROLL	2"W x 12"L x 2"H	#925822
CHEST ROLL	3"W x 6"L x 3"H	#925823
CHEST ROLL	3"W x 8"L x 3"H	#925824
CHEST ROLL	3"W x 10"L x 3"H	#925825
CHEST ROLL	3"W x 12"L x 3"H	#925826
CHEST ROLL	4"W x 12"L x 3"H	#91-387
CHEST ROLL	4"W x 16"L x 3"H	#91-388
CHEST ROLL	4"W x 20"L x 3"H	#91-389
CHEST ROLL	6"W x 12"L x 5"H	#91-390
CHEST ROLL	6"W x 16"L x 5"H	#91-391
CHEST ROLL	6"W x 20"L x 5"H	#91-392
CONTOURED CHEST ROLL	S	#925827
CONTOURED CHEST ROLL	M	#925828
CONTOURED CHEST ROLL	L	#925829
CHEST POSITIONER	LEFT	#925830
CHEST POSITIONER	RIGHT	#925831
UNIVERSAL POSITIONER	S	#95-945
UNIVERSAL POSITIONER	L	#95-946
UNIVERSAL POSITIONER	XL	#925834
DOME POSITIONER		#95-116
GRADUATED THORACIC POSITIONER		#95-113
AXILLARY ROLL COVER		#925833
CHEST ROLL COVER		#925832
LOWER EXTREMITY		
HEEL CUP		#91-386
HEEL PROTECTOR WRAP		#95-789
COMFORT HEEL CUPS		#931331
KNEE CRUTCH PAD	S	#95-942
KNEE CRUTCH PAD	M	#95-941
KNEE CRUTCH PAD	L	#95-559
CANDY CANE STIRRUP SLEEVES		#95-839
ANKLE STIRRUP STRAP PADS	S	#931334
ANKLE STIRRUP STRAP PADS	L	#938138
LITHOTOMY BOOT STIRRUP PAD		#9-244
LITHOTOMY BOOT STIRRUP PAD, SAM		#925816
TRAPEZOID POSITIONER SUPPORT		#925839

DESCRIPTION	SIZE	ITEM #
UPPER EXTREMITY		
ARMBOARD PAD	S	#930238
ARMBOARD PAD	L	#930239
ARMBOARD PAD	XL	#931333
CONTOURED ARMBOARD PAD	S	#936073
CONTOURED ARMBOARD PAD	M	#931329
CONTOURED ARMBOARD PAD	L	#936072
ULNAR NERVE PROTECTOR		#95-787
ULNAR FOREARM PROTECTOR		#95-788
HAND/WRIST SUPPORT		#936721
SHOULDER SUPPORT PADS		#925812
TOBOGGAN PAD		#935098
MAYO STAND PAD		#925793
SURFACE OVERLAYS		
VAC-PAC COVER	UNIVERSAL	#95-548
VAC-PAC COVER	T-SHAPED	#925835
FULL-LENGTH		#95-930
MEDIUM-LENGTH		#95-931
SHORT		#95-939
HOURGLASS		#925792
MEDIUM w/PERINEAL CUTOUT		#925790
THREE-PIECE FULL SET		#95-363
HEAD SECTION		#95-364
TORSO SECTION w/PERINEAL CUTOUT		#95-365
FOOT SECTION w/PERINEAL CUTOUT		#95-366
PEDIATRIC		
ISOLETTE/BASSINET PAD	12"W	#925794
ISOLETTE/BASSINET PAD	16"W	#925796
PEDIATRIC POSITIONING SET		#925795
FRAME PADS		
WILSON PAD SET	½"W	#95-837
WILSON PAD SET	1"W	#95-838
KAMBIN FRAME PAD		#925836
MONTREAL FRAME HIP PAD		#927431
MCGUIRE FRAME	ILIAC CREST PAD	#925837
MCGUIRE FRAME	FRAME PAD COVER	#925818
MCGUIRE FRAME	SACRAL PAD	#925819
RELTON FRAME PAD SET		#95-114
ANDREWS FRAME PAD	ILIAC STRAP	#925798
ANDREWS FRAME PAD	KNEE PAD	#925797
FRACTURE TABLE	HIP PAD	#925817
FRACTURE TABLE	POST PAD	#925838
PATIENT POSITIONING SETS		
BASIC POSITIONING SET		#937918
SELECT POSITIONING SET		#937917
PRONE POSITIONING SET		#937919
SUPINE POSITIONING SET		#937916
CARDIOVASCULAR POSITIONING SET		#937920



Any order placed before 3 p.m. EST, if in stock, ships the same day!
Prices, product availability, and specs are subject to change.

Pub. #7492REV0618
©2016 AliMed® Inc. All rights reserved. • Printed in USA