

# GloboPed® Forefoot Relief Shoe

## OFF-LOADING SHOE.

### Indications

- Diabetes mellitus ulcers
- Neuropathic forefoot ulcers
- Injuries of the forefoot
- Rheumatic feet

Can be worn on right and left foot.

**Three straps**  
for optimal fit and stability

**Seamless soft inner lining**

**Forefoot Shield**  
prevents undesired contact with objects

Reduced convalescent time, allowing for early patient mobility



U.S. Shoe Size	small	medium	large	extra large
Men	3½ - 5½	6 - 9½	10 - 12½	13 - 15
Women	6 - 8½	9 - 11½	12 - 13	13 +
Order Number	36627100010000	36627100020000	36627100030000	36627100040000

Sold individually