

ViscoHeel® K

VISCOELASTIC HEEL CUSHIONS FOR VARUS/VALGUS ALIGNMENT.

Indications

ViscoHeel® K

- Varus or valgus alignment
- Pronation and supination of the heel

Raised wedge can be placed medially or laterally to correct pronation/supination of the heel, or varus/valgus alignment



Viscoelastic, skin friendly material reduces shock loads and discomfort to the ankle, knee, hip, and spine

U.S. Shoe Size	1	2	3	Order no. ViscoHeel K
Men	children	3½ - 7½	8 - 12½	66682*
Women	children	5 - 9	9½ - 14	

*Specify size, Sold as a pair