




AliMed[®]
Helping people...help people™

Instructions for measuring AliMed Radiation Attenuation Gloves *Free of powder, latex*, and lead*

Measure the width of the palm of your dominant hand in millimeters or use our easy measurement guide on reverse



Sizing and Ordering Information, 5 pr/bx

Glove Size	AliMed Original		AliGuard	
	Palm Width (mm)	Item #	Palm Width (mm)	Item #
6	—	—	77-88	#960842
6½	77-88	#95-493	83-95	#960843
7	83-95	#95-494	89-101	#960844
7½	89-101	#95-495	96-107	#960845
8	96-107	#95-496	102-114	#960846
8½	102-114	#95-497	108-120	#960847
9	>114	#95-498	115-126	#960848
9½	—	—	121-132	#960849

*Not made with natural rubber latex.