

AliMed ONE[®]

Wireless Fall Monitoring System



The fully scalable **AliMed ONE® Wireless Fall Monitoring System** operates from one central Alarm Base and is easily configured to manage up to six compatible wireless monitoring devices, plus a Hall Alert.

Intuitive Operation

- Easy to learn, easy to use, easy to manage
- Unique EZPair™ system on all components

Quicker Response for Patient Focus


- Hall Alert provides visual and audible alerts
- Reset alarm system from hall to quickly eliminate in-room noise and patient stress

Easy to Clean

- Sealed Sensor Pad Transmitter is water-resistant (IPX6) with waterproof connector
- Fully gasketed speakers on Alarm Base and Hall Alert allow for complete cleaning


Wireless at a Wired Price

- Same sensor pad design as a wired system with a reusable wireless transmitter
 - Long-lasting battery life*: Alarm Base/Hall Alert, 4-month life; Transmitter, 2-year life
- *Based on normal usage



Alarm Base


Manages up to six wireless devices simultaneously, plus a Hall Alert. Integral nurse call cable with mounting bracket included.



Sensor Pad Transmitter

Connects to any compatible sensor pad (see below) for “EZ” wireless pairing with the Alarm Base. Water-resistant.


One transmitter required per sensor pad. Pairs directly with Alarm Base.



Hall Alert

Provides visual and audible hallway notifications for reduced patient stress—reset Alarm Base directly from the Hall Alert. Mounting bracket included.

Pairs directly with Alarm Base.



Motion Sensor

Adds another line of defense against elopement when motion is detected. Easily mounts to wall for protection anywhere in the room.

Pairs directly with Alarm Base.



Floor Mat

Warns staff when wandering patients attempt to leave their bed. Features a slip-resistant bottom with low-profile edges to minimize tripping.


Included transmitter pairs directly with Alarm Base.



Sensor Pads

Single-patient use, pressure-sensitive, and effective for dependable all-day monitoring in a **Bed** or **Chair**, or for security and privacy during **Toileting**.

Each sensor pad requires a transmitter.



Coming soon!

Sensor Pad Transmitter

Sensor Pad Transmitter

Sensor Pad Transmitter

Clinical, Operational, and Financial Benefits

Clinical

- Simple, easy-to-configure user interface and components
- Reduces in-room alarm noise and patient stress
- Wireless design eliminates trip hazards

Operational

- Faster response time with hallway notifications
- Allows various room configurations
- EZPair™ system for easy setup

Financial

- Easy-to-use wireless system at a wired price
- Long battery life for reliable monitoring
- Interchangeable, reusable sensor pad transmitters with a two-year life span

System Components



Alarm Base

Central unit pairs with up to 6 compatible monitoring devices • Integral 1-ft. nurse call cable for connection with nurse call system • Mounting bracket included • Requires 3 AA batteries (not included)

#713815



Hall Alert

Pairs with Alarm Base • Resets Alarm Base from unit • Mounting bracket included • Requires 3 AA batteries (not included)

#713816



Motion Sensor

Pairs with Alarm Base • Mounting bracket included • Requires one 9V battery (not included)


#713817



Sensor Pad Transmitter

Pairs with Alarm Base to connect compatible sensor pads • One required for each sensor pad • Sealed, non-replaceable battery included

#713818



Floor Mat

23¾"W x 35¾"L • Transmitter included • No-slip textured top with nonskid bottom • Tan

Contact your local AliMed representative for details.

Bed Sensor Pad



Sensor Pads

Transmitter required for each sensor pad • Single-patient use • **Bed Sensor Pad:** 10"W x 30"L • Foldable for storage • 30-day warranty • **Chair Sensor Pad:** 10"W x 15"L • 30-day warranty • **Toilet Sensor Pad:** 2"W x 5½"L • Adhesive-backed with two adhesive cord clips • 14-day warranty

#713819* Bed Sensor Pad

#713821* Chair Sensor Pad

#713823 Toilet Sensor Pad

*6-Month and 1-Year warranty options available.

Chair Sensor Pad



Toilet Sensor Pad

