


# Innovative, premium surgical positioning products



AliMed is partnering with KYRA Medical, a manufacturer of innovative surgical positioning products, to provide the next generation of lithotomy stirrups, disposables, armboards, spine positioners, and more focused on infection control and safety for O.R. staff and patients. Plus, all KYRA products are **ready to ship** so you're ready for any procedure!



- **Lithotomy**
- **Spine**
- **Orthopedics**
- **Armboards**
- **Sockets**
- **Straps and more!**

**22° MORE**  
*range of motion*  
compared to traditional stirrups  
See inside for details

## KYRA Curve™ Spine Frame Systems

- Provides a stable platform and maximum lordosis for spine surgeries
- Intraoperatively adjustable for optimal site access
- Gel pads sized to accommodate larger patients
- Easy-to-clean fabric-free silicone Secure-Lok™ straps
- Fits “Jackson style” rails and O.R. tables
- Radiolucent
- 600-lb. capacity
- System includes: Curve spine frame, gel pads, crank handle, strap connector, and storage cart

#2970021381 CURVE SPINE FRAME SYSTEM SET

Accessories and replacement parts at [AliMed.com](http://AliMed.com).

(MD) • 65 California residents, see below.

### Compare to



Allen® Bow Frame  
#A-70800

CAPACITY CAPACITY  
**600 lbs.**  
CAPACITY




Storage Cart  
for secure transport



Integrated gel pads  
allow for more  
abdomen space

## KYRA Arthroscopic Leg Holders



- Properly positions leg up to 35" in circumference
- Includes disposable pad for skin protection

#2970018050 ARTHROSCOPIC LEG HOLDER EACH

#2970018046 REPL. DISPOSABLE PAD 12/PACK

(MD) • 65 California residents, see below.

### Compare to



Allen® Arthroscopic  
Leg Holder  
#A-10015-A2

## KYRA Lateral Positioning Sets



3-Piece Set  
includes 2  
Rectangular  
and 1 Circular  
Positioner

- Secure, comfortable, and precise lateral positioning
- Swivel heads conform to patient anatomy


#2970018047 3-PIECE SET SET

#2970018049 RECTANGULAR POSITIONER EACH

#2970018048 CIRCULAR POSITIONER EACH

(MD) • 65 California residents, see below.

### Compare to



Allen® Lateral Brace  
#O-UPBL3S-A2

## KYRA Goepel Knee Crutches



- Easy to adjust for lithotomy positioning
- Pads included

#2970018034 GOEPEL KNEE CRUTCH PAIR

(MD) • 65 California residents, see below.

### Compare to



Allen® Simplicity™  
Knee Crutch  
#O-KCP

## KYRA Blue 400 Lithotomy Stirrups



CAPACITY CAPACITY  
**400 lbs.**  
CAPACITY

Fabric-free silicone  
Secure-Lok™ straps  
improve infection control

- Non-finned boot
- Same features as KYRA Comfort™ Stirrups

#2970018052 BLUE 400 LITHOTOMY STIRRUP PAIR

(MD) • 65 California residents, see below.

### Compare to




Allen® PAL Pro  
Lithotomy Stirrup  
#A-10052

California Proposition 65 Statement

65 ⚠️ WARNING: This product can expose you to chemicals including Di(2-Ethylhexyl) Phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <https://www.p65warnings.ca.gov>.

# KYRA Comfort™ Family of Lithotomy Stirrups

- Curved rod provides more range of motion and working space
- Eliminates hook-and-loop straps for improved infection control
- Larger boot and more padding supports a range of patient sizes
- Available in 350-lb., 500-lb., or 800-lb. capacities



**22° MORE**  
range of motion  
compared to traditional stirrups

Fabric-free silicone Secure-Lok™ straps improve infection control



Curved rod eliminates pinch hazards

Scan to watch the videos!



## Still using Allen Boots?

Try our Stirrup Replacement Pads at a fraction of the cost!



Models available for Yellowfins®, Ultrafins®, and non-finned stirrups

Scan to shop Stirrup Replacement Pads



CAPACITY CAPACITY  
**350 lbs.**  
CAPACITY CAPACITY

Comfort™ 350 w/Secure-Lok™ straps

#2970017788 w/SECURE-LOK PAIR  
#2970017787 w/HOOK-AND-LOOP PAIR  
MD • 65 California residents, see left page.

Compare to

Allen® Yellofins® Lithotomy Stirrup #O-YFASI



CAPACITY CAPACITY  
**500 lbs.**  
CAPACITY CAPACITY

Comfort™ 500 w/Secure-Lok™ straps

#2970017786 w/SECURE-LOK PAIR  
#2970017785 w/HOOK-AND-LOOP PAIR  
MD • 65 California residents, see left page.

Compare to

Allen® Yellofins® Elite Lithotomy Stirrup #O-YFES



CAPACITY CAPACITY  
**800 lbs.**  
CAPACITY CAPACITY

Comfort™ 800 w/Secure-Lok™ straps

#2970018053 w/SECURE-LOK PAIR  
MD • 65 California residents, see left page.

Compare to

Allen® Ultrafins® Lithotomy Stirrup #O-UFAS



## KYRA Round Post Clamps



FLAT OR ROUND

- Mounts and locks anywhere on the rail
- Autoclavable

#2970018031 ROUND POST CLAMP EACH  
MD • 65 California residents, see left page.

Compare to

Allen® Easy Lock® Socket #A-40022-A2



## KYRA Secure-Release™ Clamps



FLAT

- Prevents false rail engagement
- Mounts and locks in one movement

#2970017790 SECURE-RELEASE CLAMP EACH  
MD • 65 California residents, see left page.

Compare to

Allen® Easy Lock® Blade Clamp #A-40040-A2



# Shop these KYRA accessories...



Clean™ Lithotomy Stirrups



250 Junior Lithotomy Stirrups



Simple Lithotomy Stirrups



Lithotomy Stirrup Carts



Well Leg Holders



Lift Assist Knee Crutches



Flexible Anesthesia Screens



Secondary Shoulder Supports



Swivel Armboards



Trigger Armboards



Clamshell Armboard Pads




Table Width Extenders



Urology Attachment Bands and Drainage Bags



Patient Safety Straps



Nissen Straps



Simple Blade Clamps

...and so much more at [AliMed.com/kyra!](https://www.ali-med.com/kyra)



California Proposition 65 Statement

**65** ⚠️ **WARNING:** This product can expose you to chemicals including Di(2-Ethylhexyl) Phthalate (DEHP), which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to <https://www.p65warnings.ca.gov>.

Pub. #8092 • Product #11628

©2023 AliMed® Inc. All rights reserved.