

## Easier to clean and disinfect



- Designed to reduce cross-contamination risks

- Pivots for bilateral use

- Removable pads for easy cleaning

### Arthroscopic Leg Holder

A proven design made better by adding measures to reduce cross-contamination hazards. Removable urethane pads for ease of cleaning, with a frame constructed of stainless steel and aluminum to withstand repeated cleaning. Keeps the operative leg in place during surgical procedures—pivots 180°, for use on left or right leg. Firm pads provide stability for the thigh. Clark Socket required (FREEDOM<sup>®</sup> Clark Socket recommended, sold separately below). Adjusts 9" horizontally and 10" vertically. 12"L stainless steel mounting post. 100-lb. capacity.

### Ordering Information

#937824 Arthroscopic Leg Holder  
#934725 FREEDOM<sup>®</sup> Clark Socket



See reverse for  
additional products 

# O.R. Table Accessories



**#931962** Arthroscopic Single-Use Well-Leg Holder and Knee Kit



**#931939** Arthroscopic Reusable Vinyl-Coated Well-Leg Holder




**#934396** TKR Footrest



**#934725** FREEDOM® Clark Socket Table Clamp



**#937932** FREEDOM® 360 Socket Table Clamp



**#937344** Clark Socket Plus EZ-On Table Clamp

## Partners in the Continuum of Care

AliMed is the only niche manufacturer and distributor with on-site customization capabilities. This complements a vast portfolio of in-stock, specialty healthcare products designed to improve patient outcomes and experiences. Serving patient needs spanning hospital to home, AliMed has built its reputation on responsiveness, exceptional customer service, and high-quality offerings.

### The AliMed Advantage:

- Field sales and dedicated customer service support
- Same-day shipping for in-stock items\*
- Efficient ordering through EDI
- E-commerce site and specialized catalogs to easily find and order products
- On-site customization capabilities
- Competitive pricing through contracts

\*Place your order before 3 p.m. EST—if it's in stock, it ships the same day!  
Prices, product availability, and specs are subject to change.

Pub. #7173REV0318 Product #11139  
©2018 AliMed® Inc. All rights reserved. • Printed in USA.